UNCLASSIFIED

AD NUMBER AD468375 **NEW LIMITATION CHANGE** TO Approved for public release, distribution unlimited **FROM** Distribution authorized to U.S. Gov't. agencies and their contractors; Administrative/Opertional Use; May 1965. Other requests shall be referred to USAF School of Aerospace Medicine, Brooks AFB, TX. **AUTHORITY** SAM ltr, 8 Jun 1966

SECURITY MARKING

The classified or limited status of this report applies to each page, unless otherwise marked.

Separate page printouts MUST be marked accordingly.

THIS DOCUMENT CONTAINS INFORMATION AFFECTING THE NATIONAL DEFENSE OF THE UNITED STATES WITHIN THE MEANING OF THE ESPIONAGE LAWS, TITLE 18, U.S.C., SECTIONS 793 AND 794. THE TRANSMISSION OR THE REVELATION OF ITS CONTENTS IN ANY MANNER TO AN UNAUTHORIZED PERSON IS PROHIBITED BY LAW.

NOTICE: When government or other drawings, specifications or other data are used for any purpose other than in connection with a definitely related government procurement operation, the U. S. Government thereby incurs no responsibility, nor any obligation whatsoever; and the fact that the Government may have formulated, furnished, or in any way supplied the said drawings, specifications, or other data is not to be regarded by implication or otherwise as in any manner licensing the holder or any other person or corporation, or conveying any rights or permission to manufacture, use or sell any patented invention that may in any way be related thereto.

PJ 888 375-

TELUS (Telemetric Universal Sensor)

WILLIAM G. GLENN, Ph.D.
WESLEY E. PRATHER
HEINZ A. JAEGER

FOR FEDE TECHNI	ARINGHO RAL SCLEN ICAL INFOR	TIFIC AND
Bardcopy	Microfiche	,
\$1.60	\$,50	14 pp 78
ARC	HIVE C	OPY

May 1965

PROCESSING COPY

USAF School of Aerospace Medicine Aerospace Medical Division (AFSC) Brooks Air Force Base, Texas

TELUS

(Telemetric Universal Sensor)

WILLIAM G. GLENN, Ph.D.
WESLEY E. PRATHER
HEINZ A. JAEGER

FOREWORD

This report was prepared in the Microbiology Section of the Biosciences Branch under Task No. 775402. The paper was submitted for publication on 18 December 1964. The work was accomplished between July 1963 and December 1964.

The coder and searcher used in this study were purchased with funds from the Laboratory Director's Fund of the USAF School of Aerospace Medicine.

This technical documentary report has been reviewed and is approved.

Harold V. Ellingson HAROLD V. ELLINGSON

Colonel, USAF, MC

Commander

ABSTRACT

There was need for a flexible complex of receiving and evaluating instruments for sensing physiologic and biologic analyses performed remotely by various field transducers. This need has been met by the design and development of TELUS (Telemetric Universal Sensor), a one-man laboratory console capable of receiving, quantitating, comparing, coding, storing, searching, retrieving, and distributing electrical and electromagnetic data. These data are received and distributed by radio and telephone. Field tests of TELUS indicate good performance characteristics for evaluating four channels of telemetered data and providing two-way communication between the laboratory and remote testing areas.

TELUS

(Telemetric Universal Sensor)

I. INTRODUCTION

When the feasibility of short- and longrange telemetry from certain analytic field sensors for biologic fluids was established (1-6), it became evident that an integrated functional console to be operated in the laboratory was required. The console was to receive, quantitate, compare, code, store, search, retrieve, and distribute data from field sensors located at some remote site. These functional needs generated the design and development of the TELUS (Telemetric Universal Sensor). TELUS provides the capability for monitoring any electrical or electromagnetic output from field-sensing instruments and quantitating these data into readout instruments (oscilloscope, strip-chart recorder), while simultaneously distributing the data signals to other laboratories.

II. SUMMARY

Successful integration of modified, commercially obtainable components has resulted in the functional design and development of a laboratory console that is suitable for sensing several physiologic and biologic transducers located remotely. Field tests of TELUS have resulted in a telephone-radio, two-way communication system that can be operated by one man in a laboratory. The four-channel input of TELUS can be in electrical or electromagnetic form and the output distributed, taped, and viewed simultaneously on a dual-channel oscilloscope and a strip-chart recorder. Comparison of data is anticipated by use of two instrumentation tape recorders simultaneously and two-way radio telemetry from field sensors located remotely.

III. DESIGN AND DEVELOPMENT

TELUS was designed by the authors and developed entirely at the School of Aerospace Medicine from commercial components. Some of the instruments incorporated into the console were modified electronically to meet the required performance integration. The wooden console is 2 meters high, 1.5 meters wide, and 0.9 meter deep, with instruments grouped into three vertical functional modules, each of which is in a double oblique plane for easy one-man operation. Access to the back of the instruments is through three separate panel openings in the rear of the console. Two instrumentation tape recorder-reproducers flank the console in matching two-tone wooden cabinets.

IV. FUNCTIONAL CAPABILITIES

TELUS is designed to:

- a. Receive telemetered signals directly or indirectly from air or ground by radio, telephone, and magnetic tape.
- b. Initiate or relay data to space or ground stations.
 - c. Decode frequency-modulated signals.
- d. Place a binary code (000-999) on both incoming and outgoing data and utilize this code through an automatic searching system to retrieve information from magnetic tape.
- e. Present channels of data for viewing on a dual-channel oscilloscope and a two-channel strip-chart recorder.

f. Simultaneously receive four channels of the same or different information (either frequency-modulated or direct—i. e., voice). This information can be simultaneously recorded on any of fourteen channels on magnetic tape.

V. APPLICATIONS

TELUS can be linked directly or indirectly to sensors for quantitation of:

- a. Physiologic data (e.g., electrocardiography, electroencephalography, respirometry, sphygmomanometry).
- b. Physical data (e.g., gas, pressure, strain, temperature).
- c. Physical-chemical data (e.g., pH, turbidimetry, spectrophotometry, reflectometry, nephelometry, Liesegang phenomenon).

VI. GENERAL OPERATION

High-frequency electromagnetic waves picked up by the antenna are conducted to the receiver (B of fig. 1). Audio signals can be received via the telephone circuit (T). Whenever necessary, the output of the receiver (or telephone) may be increased by connecting the receiver (or telephone) to one of the amplifiers (C) before the information is relayed to the next circuit. The signals are decoded by discriminators (A), according to the frequencies being used, and transposed into individual characteristics with analog or digital outputs, which can be visualized on a dual-channel oscilloscope (I) and a two-channel, galvanometric strip-chart recorder (S). Simultaneously, the signals can be stored on one to fourteen channels of magnetic tape (X), paralleled with a preselected binary code (M) and audio recording from headsets (J, J_1) . A high-frequency counter (R) can be used to visually indicate the frequency of the incoming or outgoing audio information before distribution. When the binary code used to identify data is fed onto magnetic tape, such code can be electronically sensed on playback by the tape searcher (L), which automatically stops the recorder when the desired identifying number appears in the reading head of the recorder (X). Controls (E-H) provide remote operation of one tape recorder from the right side of TELUS.

TELUS is capable of receiving frequencies of 50 to 260 mc. These signals can have bandwidths of 0 to 300 kc. Received signals can be of the video, audio, or subcarrier type (400 cps to 80 kc./sec.). Radio-frequency carriers can be amplitude- or frequency-modulated.

Distribution of data through and within TELUS is accomplished by interconnecting components in various sequences (modes). Central programing and distribution were selected to maximize the utilization of available instruments. These essential instruments are assembled in the console with outputs and inputs connected to a central grid system. This central grid-distribution system is adjusted for 1.4-volt distribution into 1,000-ohm loads to allow multiple bridging for chart, oscilloscope, recorder, and amplifier circuits, while minimizing cross-talk and loading problems. TELUS is provided with forty vertical grids and thirty horizontal grids (O of fig. 1). The input or output of each instrument is assigned a position at one of the forty vertical grids. The thirty horizontal grids are called switching channels. Any of the vertical instrument connections of TELUS can be connected to any of the horizontal switching channels by placing a peg in the grid where the circuits cross. The horizontal switching circuits are brought out to the rear terminals of TELUS. Antenna, telephone lines, and tape recorders are connected to these rear terminals. For each specific experiment there is a corresponding arrangement of pegs in the grid-distribution block. A plastic overlay is made with holes where the pegs are to be inserted. To preset TELUS for a given mode of operation, the proper overlay is selected and pegs inserted. This gives TELUS the flexibility of a universal system with the reliability of an integrated system.

TELUS is organized around a system of FM/FM telemetry. An advantage of this system is that data can be provided for a new




FIGURE 1

TELUS (Telemetric Universal Sensor).

- A. Discriminator (decoding) rack
- B. Receiver (50 to 260 mc.)
- C. Operational amplifiers
- D. Filters
- E. Tape recorder (remote control)
- F. Headset amplifier
- G. Delayed playback control
- H -- 'I1. Headset interconnect
- I. Oscilloscope (dual-channel)
- $J = J_1$. Headset
- K. Monitor meter
- L. Magnetic tape searcher (digital)

- M. Magnetic tape coder (digital)
- N. Switch-pin rack
- O. Program boards (switch-pin type)
- P. Test panel (terminal pin jacks under the cover)
- Q. Test panel selector switches
- R. Audio frequency counter
- S. Strip-chart recorder (dual-channel)
- T. Telephone (send, receive)
- U. Internal timer
- V. Telephone selector switches
- W. Telephone key switch
- X. Magnetic tape recorder (two separate recorders, fourteen channels each, connected into TELUS)

measurement by the simple addition of another subcarrier, Data come from instruments widely separated; thus, the subcarrier system is the best method of combining data.

TELUS includes a transmitter that can relay any signal via other communication links over 245.8 mc. FM (approx. 10 watts).

VII. GRAPHIC INTERPRETATION

When the feasibility of transmitting data from a prototype Bio-telescanner via the Military Affiliated Radio Service from one hemisphere to another was demonstrated in 1962 (3, 6), emphasis was placed on refinement and integration of the Bio-courier (2)-TELUS combination.

One type of biologic analyses, precipitin (antigen-antibody) systems in agar gels, is shown in figure 2. Columns of this type are photoelectrically scanned by the Bio-telescanner (5) shown in figure 3. Changes in the photoelectric cell output are fed to the Bio-transceiver (fig. 4), where the transmitter is located. The transceiver was designed as a companion instrument to provide sustained battery power for the telescanner and a two-way telemetry link between an operator in the field and the scientist in the laboratory. (For space vehicles, a transceiver would not be required since the 28-volt input to the scanner and the telemetry would be supplied by the vehicle)

Quantitation of precipitin reactions and other essential data sent from the Bio-telescanner in the field to TELUS takes the graphic form shown in figure 5. In this instance, data from two subcarriers (1.7 kc. and 1.3 kc.) are shown. The upper half of figure 5 (1.7 kc.) shows a series of three square waves on the The amplitude of these indicates the temperature in the telescanner and the first and second digits of the reaction being scanned. Should atmospheric disturbances (static) affect transmission of the data, these would be shown as erratic spikes in the horizontal components of the square waves. An indicator of this nature is essential for determining the quality of the transmission. The sequence of spikes further to the right on the upper half of figure 5 represents the binary (yes-no) code placed on each scan by the TELUS operator. The 12 spikes are divided into three groups of 4. Within each group the spikes are given the value of 1-2-4-8, reading from right to left. Where spikes are toward the top (position yes-e.g., 2, 1), the values for each group are added to give the total for that group. By combinations of yes and no within each group, a three-digit numeric code from 000 to 999 can be registered. Binary codes are placed on the 1.7-kc. signals by the TELUS operator, thus, energizing the magnetic tape coder (M) of fig. 1) to be used later by the magnetic tape searcher (L of fig. 1) for automatic searching and retrieval.

The lower half of figure 5 is the 1.3-kc. subcarrier tracing the graphic picture of a reaction column. An explanation of the graph is given in the figure legend.

VIII. FIELD TESTS

The console described in this report is used routinely to monitor precipitin analyses via telemetry from an automatically programed Bio-telescanner (5). In addition, during the past year, TELUS has been used several times in a two-way communication system to monitor physiologic changes of subjects in altitude chambers. In these instances, communication from TELUS was by telephone and data from the chambers sent via radio to the laboratory. Other field tests, as exercises, have been done in two-way communication by a combination of radio and telephone linking the Bio-telescanner and its operator, located in an altitude chamber, with TELUS located in another building about 200 yards away.

Successful operation of the Bio-courier—TELUS telemetry system will depend on established communication links. Anticipated paths and field areas to be utilized are shown in figure 6. Each of the remote locations has numerous microbiologic problems that can be solved better and more efficiently through cooperative research via telemetry. It is this approach for which TELUS was designed.

FIGURE 2

Precipitin (antigen-antibody) reactions in agar diffusion columns. The glass tubes are 60 mm. long and 4 mm., outside diameter. In most instances, photometric quantitation of the white zones at an early time when they are aggregated close to the interface gives the essential information.



FIGURE 3

Bio-telescanner, model X-63B. Reaction columns are positioned in the round, black magazine. The black tube to the right moves the scanning beam in a vertical direction. In the later model, X-63C, the electronic deck in the foreground is positioned in a false bottom under the scanner.





FIGURE 4

Bio-telescanner, model X-63C (left) and Bio-transceiver model X-63BR. The transceiver provides sustained battery power to the scanner and permits two-way communication between the operator in the field and the scientist at TELUS in the laboratory.







FIGURE 5

Graphic form of the telemetry data received by TELUS from the Bio-telescanner. The upper half of the figure (1.7 kc.) shows the temperature, magazine position, and an identifying binary code as explained in the text. The lower half (1.3 kc.) represents the following information concerning the reaction columns being scanned: A, magazine changing position; B, deflection where the antigen (see fig. 2) meets the air; C, reference mark on the column; D and H, interface where antigen meets antiserum-agar; E and G, precipitin zone; F, antiserum-agar without reactions.



M.A.S.S. Base Station

FIGURE 6

Anticipated sites for field studies and the communication links in the Biocourier-TELUS system of telemetry of biologic reactions.

REFERENCES

- Glenn, W. G. Comparative immunobiology takes to the air. Serological Museum Bulletin No. 28. New Brunswick, N. J.: Rutgers University, Nov. 1962.
- 2. Glenn, W. G. The Bio-courier. Comparative immunobiological support for planetary, space, and interplanetary explorations. Lectures in Aerospace Medicine, USAF School of Aerospace Medicine, Feb. 1963.
- 3. Glenn, W. G., and W. T. Prather. Immunobiological support for manned orbital and interplanetary explorations. Aerospace Med. 34:408 (1963).
- 4. Glenn, W. G., H. A. Jaeger, and W. E. Prather. Bio-instrumentation and telemetry for immuno-

- chemical analysis. Proceedings, 1963 National Telemetering Conference, Albuquerque, N. Mex. 20 May 1963.
- Glenn, W. G., H. A. Jaeger, and W. E. Prather.
 The Bio-telescanner: An instrument for telemetric quantitation of immunodiffusion (antigen-antibody) reactions. Proceedings of 1st National ISA Biomedical Sciences Instrumentation Symposium, Los Angeles, Calif., 18 June 1963.
- 6. Glenn, W. G. East meets West. The Bio-courier project. Phase I completed. Serological Museum Bulletin No. 29. New Brunswick, N.J.: Rutgers University, June 1963.

Security Classification

DOCUMENT COMPANDED TO THE PROPERTY OF THE PROP	NTROL DATA - R&D	red when th	ne overall report is classified)		
USAF School of Aerospace Medicine		Za REPORT SECURITY CLASSIFICATION Unclassified			
Aerospace Medical Division (AFSC)		26 GROUP			
Brooks Air Force Base, Texas					
3 REPORT TITLE					
TELUS (Telemetric Universal Sensor)					
4. DESCRIPTIVE NOTES (Type of report and inclusive dates) Initial report July 63 - Dec. 64					
S AUTHOR(S) (Last name, first name, initial) Glenn, William G.					
Prather, Wesley E.					
Jaeger, Heinz A.					
6: REPORT DATE May 65	7ª TOTAL NO OF PAG	G F. S	75 NO OF REFS		
SE CONTRACT OR GRANT NO.	90 ORIGINATOR'S REPORT NUMBER(S)				
b. PROJECT NO	SAM-TR-65-1				
Task No. 775402					
	95 OTHER REPORT NO this report)	U(3) (Any	offier numbers that may be easigned		
4					
10 AVAILABILITY/LIMITATION NOTICES					
Qualified requesters may obtain copies	of this report	from DI	DC.		
11 SUPPLEMENTARY NOTES	12 SPONSORING MILITA	ARY ACTIV	VITY		
		_	space Medicine		
	Aerospace Medical Division (AFSC)				
	Brooks Air Fo				
13 ABSTRACT					

There was need for a flexible complex of receiving and evaluating instruments for sensing physiologic and biologic analyses performed remotely by various field transducers. This need has been met by the design and development of TELUS (Telemetric Universal Sensor), a one-man laboratory console capable of receiving, quantitating, comparing, coding, storing, searching, retrieving, and distributing electrical and electromagnetic data. These data are received and distributed by radio and telephone. Field tests of TELUS indicate good performance characteristics for evaluating four channels of telemetered data and providing two-way communication between the laboratory and remote testing areas.

DD 5084. 1473

Unclassified

Security Classification 14. KEY WORDS	LIN	KA	LINK B		LINKC	
	ROLE	WT	ROLE	WT	ROLE	WT
Instrumentation						
Instrumentation						
Telemetry						,
Biologic sensing						

INSTRUCTIONS

- ORIGINATING ACTIVITY: Enter the name and address of the contractor, subcontractor, grantee, Department of Defense activity or other organization (corporate author) issuing the report.
- 2a. REPORT SECURITY CLASSIFICATION: Enter the overall security classification of the report. Indicate whether "Restricted Data" is included. Marking is to be in accordance with appropriate security regulations.
- 2b. GROUP: Automatic downgrading is specified in DoD Directive 5200.10 and Armed Forces Industrial Manual. Enter the group number. Also, when applicable, show that optional markings have been used for Group 3 and Group 4 as authorized.
- 3. REPORT TITLE: Enter the complete report title in all capital letters. Titles in all cases should be unclassified. If a meaningful title cannot be selected without classification, show title classification in all capitals in parenthesis immediately following the title.
- 4. DESCRIPTIVE NOTES: If appropriate, enter the type of report, e.g., interim, progress, summary, annual, or final. Give the inclusive dates when a specific reporting period is covered.
- 5. AUTHOR(S): Enter the name(s) of author(s) as shown on or in the report. Enter tast name, first name, middle initial. If military, show rank and branch of service. The name of the principal author is an absolute minimum requirement.
- 6. REPORT DATE: Enter the date of the report as day, month, year, or month, year. If more than one date appears on the report, use date of publication.
- 7a. TOTAL NUMBER OF PAGES: The total page count should follow normal pagination procedures, i.e., enter the number of pages containing information.
- 7b. NUMBER OF REFERENCES: Enter the total number of references cited in the report.
- 8a. CONTRACT OR GRANT NUMBER: If appropriate, enter the applicable number of the contract or grant under which the report was writter.
- 8b, 8c, & 8d. PROJECT NUMBER: Enter the appropriate military department identification, such as project number, subproject number, system numbers, task number, etc.
- 9a. ORIGINATOR'S REPORT NUMBER(S): Enter the official report number by which the document will be identified and controlled by the originating activity. This number must be unique to this report.
- 9b. OTHER REPORT NUMBER(S): If the report has been assigned any other report numbers (either by the originator or by the aponsor), also enter this number(s).
- 10. AVAILABILITY/LIMITATION NOTICES: Enter any limitations on further dissemination of the report, other than those

imposed by security classification, using standard statements such as:

- (1) "Qualified requesters may obtain copies of this report from DDC."
- (2) "Foreign announcement and dissemination of this report by DDC is not authorized."
- (3) "U. S. Government agencies may obtain copies of this report directly from DDC. Other qualified DDC users shall request through
- (4) "U. S. military agencies may obtain copies of this report directly from DDC. Other qualified users shall request through
- (5) "All distribution of this report is controlled. Qualified DDC users shall request through

If the report has been furnished to the Office of Technical Services, Department of Commerce, for sale to the public, Indicate this fact and enter the price, if known

- 11. SUPPLEMENTARY NOTES: Use for additional explanatory notes.
- 12. SPONSORING MILITARY ACTIVITY: Enter the name of the departmental project office or laboratory sponsoring (paying for) the research and development. Include address.
- 13. ABSTRACT: Enter an abstract giving a brief and factual summary of the document indicative of the report, even though it may also appear elsewhere in the body of the technical report. If additional space is required, a continuation sheet shall be attached.

It is highly desirable that the abstract of classified reports be unclassified. Each paragraph of the abstract shall end with an indication of the military security classification of the information in the paragraph, represented as (TS), (S), (C), or (U).

There is no limitation on the length of the abstract. However, the suggested length is from 150 to 225 words.

14. KEY WORDS: Key words are technically meaningful terms or short phrases that characterize a report and may be used as index entries for cataloging the report. Key words must be selected so that no security classification is required. Identifiers, such as equipment model designation, trade name, military project code name, geographic location, may be used as key words but will be followed by an indication of technical context. The assignment of links, rules, and weights is optional.